










---






# Soffio PLUS

---

A

**ENGLISH** pag. 5

**FRANÇAIS** pag. 11

**DEUTSCH** pag. 18

**ESPAÑOL** pag. 25

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 Please read the instructions and the important notes (1.3) carefully before using the device.

## 1 - IMPORTANT INDICATIONS

### 1.1 - APPLICATION

SOFFIO PLUS is an aerosol therapy device for administering drugs in spray form.

SOFFIO PLUS features a series of innovative solutions that make it effective for treating breathing disorders.

The PERFECTA kit (nebuliser + spacer mask) and mouthpiece are effective in treating disorders of the lower airways (including bronchial asthma, chronic bronchitis, pulmonary emphysema and cystic fibrosis).

 The administration of drugs requires a medical prescription specifying the type of drug, the dosage and the length of treatment.

This aerosol therapy device and its accessories have been designed and made in Italy in accordance with EU Directive 93/42/EC regarding medical devices (and subsequent updates) and bear the CE mark.

This product complies with the restriction on the use of certain hazardous substances in electric and electronic devices (European Directive 2011/65/EU).

Device suitable for continuous operation.

### 1.2 - SYMBOLS USED

Class II apparatus

Apparatus with type BF applied part

Important - read the instructions carefully

Alternating current

Apparatus ON "I"

Apparatus OFF "O"

The apparatus must be dismantled at the end of its lifecycle and disposed of in special bins. Do not dispose of in ordinary town waste bins (EU Directive 2002/96/EC)

Device protected against finger access to hazardous parts and drops of water falling from above (IP21).

Keep the device dry

Manufacturer

Does not contain phthalates

No traces of natural latex

### 1.3 - IMPORTANT NOTES

- Read the instructions carefully and keep them in a safe place for future reference.
- Do not use the apparatus or accessories for anything other than the intended purpose.
- Children must only be allowed to use the device under the supervision of an adult who has read and understands the instructions.
- Some components are small enough to swallow and may lead to suffocation.
- The power cable and connecting pipe are long enough to pose a strangulation hazard.
- Read the list of components in the table in section 7 to see if there are any components that have caused an allergic reaction in the past.
- Do not use essential oils (e.g. with menthol or eucalyptus) as they are not compatible with the material of which the nebuliser is made.
- During the session, the apparatus must stand on a firm flat surface.
- Never cover the air vents during operation.
- Keep well away from sources of heat.
- Do not touch the plug or the electrical components wet or with damp hands.
- Device protected against finger access to hazardous parts and drops of water falling from above (IP21).
- Do not pour water or other liquids onto the device.
- Do not use in the bath or shower.
- If the apparatus should all into liquid, unplug before removing it. Do not use the apparatus. Send it to an approved service centre for an overhaul.
- Beware of compressed air coming out of the air outlet as this may cause a hazard; it must only be used to supply the nebuliser.
- The apparatus is not suitable for use in the presence of anaesthetic mixtures with air, oxygen or nitrogen protoxide.
- Performance of the apparatus is only guaranteed if all the accessories used are original ones. Only original spare parts must be used.
- The manufacturer declares that the device has a useful life of at least 5 years.

## 1.4 - DESCRIPTION OF THE APPARATUS

The SOFFIO PLUS aerosol therapy device comprises a compressor operated by an electric motor in the base, a PERFECTA kit (nebuliser and volumetric mask), a mouthpiece (which can be mounted on the nebuliser instead of the mask) and a flexible tube connecting the nebuliser and the base unit.

## 1.5 - PATIENT GROUPS

The SOFFIO PLUS aerosol therapy device provides highly effective inhalation therapy for patients of any age, from young children to the elderly.

## 1.6 - CONTENTS

The package contains (Fig. A):

The basic unit (1), a power cable (2), an ON/OFF switch (3), a fuse holder and fuse (4) integrated in the power plug (5), filter slot and filter (6), an air intake (7) and nebuliser support (8).

It also contains:

- a Perfecta kit containing the following:
  - a nebuliser (upper section (9), atomiser (10), lower section (11))
  - a spacer mask (mask body (12), fitting with discharge valve (13))
- a mouthpiece (14)
- a connecting tube (15)
- carrying case (16)
- spare air filters (17)

## 2 - PREPARATION

### 2.1 - PREPARING THE BASIC UNIT

Before use, make sure the voltage indicated on the rating plate underneath the apparatus corresponds to the mains voltage and the wiring system complies with the current regulations.

Do not use if the cable is damaged.

Check that the casing is integral and there is no sign of dents or damage.

Make sure the switch (3) is in position 0 (OFF). Plug the power cable into the socket in the apparatus, insert the electric plug in the mains socket, press position I (ON) on the switch and check that air comes out of the air intake (7). Then press position 0 (OFF) on the switch.

### 2.2 - PREPARING THE PERFECTA KIT

Please read section 3.2 before using the device.

Proceed as follows to use the Perfecta kit:

- Make sure the atomiser (10) is correctly mounted on the central sleeve in the lower section (11) of the nebuliser.
- The lower section (11) has a graduated scale that shows the approximate amount of liquid inside. The purpose is merely to give an approximate indication of the contents. Use a graduated syringe to get an accurate dose. Pour the liquid into the tank (Fig. B), keeping to the dosage specified by your doctor.

Readings must be taken on the graduated scale with the atomiser mounted on the lower section.

- Screw the upper section (9) of the nebuliser onto the base
- Fit the spacer mask (13) onto the upper section (9) of the nebuliser (Fig. C)
- Connect the lower end of the tube (15) to the nebuliser and the other end to the air intake (7)
- Press position I (ON) on the switch (3)
- During treatment, hold the nebuliser upright in your hand to enable it to work efficiently
- If the above operations have been performed correctly, a spray will come out of the mask

When used by children, the spacer mask (12) must be placed gently on the face so that it covers the mouth and nose. When used by adults, it must be placed gently over the mouth to cover it. The nose must remain outside the mask. For the treatment to be effective, you should breathe slowly and deeply. If the user breathes out through the mouth, the valve on the fitting (13) releases the exhaled air (Fig. D).

Do not obstruct the valve on the fitting during use (Fig. D).

When you have used up all the liquid or wish to interrupt treatment, press (0) (OFF) on the switch.

A small amount of liquid may remain in the nebuliser at the end of the session but this has no importance.

If a small amount of moisture remains inside the air tube at the end of the session, leave the device switched on and allow it to dry by ventilation.

The presence of liquid inside the tube is the result of water vapour forming under certain environmental conditions.

### 2.3 - USING THE MOUTHPIECE

Please read section 3.3 before using the device.

To use the mouthpiece, fit it onto the upper section of the nebuliser (Fig. C).

The mouthpiece must be held in the mouth, between the teeth with the lips closed, even when you breathe out. Breathe in through the mouth and out through the nose.

During use, do not obstruct the air outlet hole in the base.

Atomisation is very rapid and the session only lasts a few minutes.

As stated under Important Notes, some liquids are not compatible with the materials of which the nebuliser, mask and the Rinowash are made. They may interact and cause damage (the materials are shown in section 7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS).

## 3 - CLEANING, DISINFECTION AND STERILISATION

Before proceeding, unplug the device and disconnect the nebuliser, Rinowash and tube.

If it is necessary to clean, disinfect or sterilise this medical device using a method other than that given in this section, contact the manufacturer first.

### 3.1 - CLEANING, DISINFECTING AND STERILISING THE BASIC UNIT

#### 3.1.1 - CLEANING THE BASIC UNIT

The basic unit (1) must not be cleaned in water as it could get inside and damage the electrical components, which would cause a safety hazard. Clean the basic unit (1) with a soft damp cloth, with a small amount of mild (non-abrasive) detergent, if necessary.

#### 3.1.2 - DISINFECTING THE BASIC UNIT

The basic unit can be disinfected using a soft sponge or cloth and some cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions.

Make sure the solution is compatible with the material of which the basic unit is made. Refer to section 7 - DISPOSING OF THE APPARATUS AND ITS COMPONENTS.

#### 3.1.3 - STERILISING THE BASIC UNIT

The basic unit (1) must not be sterilized.

### 3.2 - CLEANING, DISINFECTING AND STERILISING THE NEBULA SPACER KIT

To prevent the risk of microbial contamination, the nebuliser and mask must be cleaned and disinfected, and if necessary sterilised, after each treatment and before use, following the instructions carefully.

Before, cleaning, disinfecting or sterilising any of the components, disassemble them as follows:

- Remove the tube (15).
- Remove the mask (12, 13) from the nebuliser.
- Detach the fitting (13) from the mask (12).
- Unscrew the upper section (9).
- Remove the l'atomizzatore (10).

#### 3.2.1 - CLEANING THE PERFECTA KIT

The disassembled nebuliser (9, 10 and 11) and the disassembled spacer mask (12 and 13) can be washed in water at an approximate temperature of 40°C with the addition of detergent (refer to the manufacturer's instructions), then rinsed in running water.

When all the parts have been cleaned, dry them with a soft cloth and place in the storage compartment.

#### 3.2.2 - DISINFECTING THE PERFECTA KIT

After washing, all the components of the nebuliser (9, 10 and 11) and the spacer mask (12, 13) can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions.

The components of the nebuliser and spacer mask can be disinfected by boiling in a saucepan of water for a maximum of 10 minutes. Make sure the components are not in direct contact with the bottom of the saucepan.

Boiling may alter the colour or appearance of the components, but this will not affect operation of the device.

#### 3.2.3 - STERILISING THE PERFECTA KIT

The components of the nebuliser (9, 10, 11) and the spacer mask (12, 13) can be sterilised in an autoclave at 121°C for at least 15 minutes (but no more than 30 minutes), for up to 20 times.

Sterilisation is recommended before each session when the device is used in a hospital environment or by more than one patient.

### 3.3 - CLEANING, DISINFECTING AND STERILISING THE MOUTHPIECE

To prevent the risk of microbial contamination, the mouthpiece (14) must be cleaned and disinfected, and if necessary sterilised, after each treatment and before use, following the instructions carefully.

#### 3.3.1 - CLEANING THE MOUTHPIECE

The mouthpiece can be washed in water at an approximate temperature of 40°C with the addition of detergent (refer to the manufacturer's instructions), then rinsed in running water.

When it has been cleaned, dry it with a soft cloth and place in the storage compartment.

### 3.3.2 - DISINFECTING THE MOUTHPIECE

After washing, the mouthpiece can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions.

### 3.3.3 - STERILISING THE MOUTHPIECE

The mouthpiece can be sterilised in an autoclave at 121°C for at least 15 minutes (but no more than 30 minutes), for up to 20 times.

### 3.4 - CLEANING, DISINFECTING AND STERILIZING THE TUBE

The tube (15) needs to be cleaned after each session. Detach it from the nebuliser and leave the apparatus on until all signs of moisture have been eliminated. If the tube is particularly dirty, it needs to be replaced.

 The tube (15) must not be sterilized hot in an autoclave or boiled.

## 4 - MAINTENANCE

The SOFFIO PLUS aerosol therapy device does not require any particular maintenance. The compressor does not need to be lubricated as the piston is the dry-operating type.

A special external filter keeps out dust and impurities, and should be replaced periodically.

 Always unplug from the mains before carrying out the following operations.

### 4.1 - REPLACING THE FUSE (Fig. E)

If the apparatus does not work even when plugged in, it may be necessary to change the fuse.

Unplug the apparatus.

Use a small screwdriver (not included) to remove the fuse holder above the power socket on the apparatus.

Remove the fuse and replace it with a new one having the same features (see Section 6).

Insert the fuse holder in the slot. Plug in and check to see if it switches on.

### 4.2 - REPLACING THE AIR FILTER (Fig. F)

The air filter needs to be replaced once a year if the device is used frequently.

The air filter (6) is located underneath the appliance. Remove the cap using a screwdriver and change the air filter.

The filter cannot be washed or regenerated. It must always be replaced with a new one.

### 4.3 - MAINTENANCE OF THE ACCESSORIES

When cleaning, disinfecting or sterilising the components, you can check them visually and determine when they need to be replaced. Repeated sterilisation may affect the materials of which the accessories are made, so the components should be inspected regularly. The accessories must be replaced when there are clear signs of deterioration.

## 5 - PROBLEMS, POSSIBLE CAUSES AND SOLUTIONS

### TROUBLE-SHOOTING CHART

Problems	Possible causes	Solutions
The apparatus does not work.	<ul style="list-style-type: none"><li>- The plug is not inserted in the socket.</li><li>- The switch (3) is in position 0 (OFF).</li><li>- The fuse has blown.</li><li>- The thermostat has cut in.</li></ul>	<ul style="list-style-type: none"><li>- Insert the plug in the socket and press position I (ON) on the switch (3).</li><li>- Press position I (ON) on the switch (3).</li><li>- Replace the fuse (see 4.1).</li><li>- Clear the air vents.</li></ul>
The nebuliser does not generate a spray.	<ul style="list-style-type: none"><li>- The apparatus is switched off.</li><li>- There is no liquid in the nebuliser.</li><li>- Too much liquid has been poured into the nebuliser.</li><li>- The nebuliser has not been assembled correctly.</li><li>- The air outlet (11) in the base is blocked.</li><li>- The air tube (15) is not connected securely.</li></ul>	<ul style="list-style-type: none"><li>- Switch on the apparatus.</li><li>- Pour liquid into the nebuliser.</li><li>- Reduce the quantity of liquid in the nebuliser to below the maximum level.</li><li>- Check that the atomiser (10) is properly inserted in the base of the nebuliser (11). Screw the upper section (9) onto the base and check that air comes out of the little hole at the top of the base.</li><li>- Unblock the hole.</li><li>- Connect the tube (15) securely to the air intake on the device (7) and the one on the nebuliser.</li></ul>
Drug does not reach the mouth or there is difficulty exhaling (problem due to the spacer mask)	<ul style="list-style-type: none"><li>- The spacer mask (12, 13) has not been assembled correctly.</li><li>- The valve (Fig. D) has jammed.</li></ul>	<ul style="list-style-type: none"><li>- Reassemble the mask (12, 13) correctly.</li><li>- Release the valve.</li></ul>


If you have failed to solve the problem after performing the above checks and the device still does not work, have it checked by an authorized service centre.

## 6 - TECHNICAL DATA

Supply voltage	See rating plate on the device
Operation	Continuous
Maximum air pressure	2,5 bar
Maximum air flow	12 l/min
Dynamic air flow with nebuliser	6 l/min
Maximum input	115 VA
Protection fuse	F1,6 A
Dimensions of the basic unit	190x130x120H mm
Weight of the basic unit	1,5 kg
Noise level at 1 m	59,3 dBA (to UNI EN 13544 -1)
Minimum nebuliser filling volume	1ml
Maximum nebuliser filling volume	8 ml
Aerosol output with nebuliser (*)	0.75 ml (**)
Minimum aerosol output rate with nebuliser (**)	0.30 ml/min (**)
MMAD with nebuliser, graph 1 (*) (***)	1.9 $\mu$ m (**)

**GRAPH 1**

The MMAD values shown are measured using a MACH 2 API AEROSIZER.

The y-axis of the graphs represent the percentage of volume of atomised solution with sizes lower than the corresponding dimensional value expressed in  $\mu$ m on the x-axis.

(\*) For a quantity of 2 ml.

(\*\*) The values shown refer to the use of physiological solution (0.9% NaCl) and may vary with the drug used. The values given do not apply to drugs in suspension or high-viscosity drugs. The manufacturer of the drug can provide the necessary details.

(\*\*\*) 50% by volume of the atomised particles are less than the MMAD (Median Mass Aerodynamic Diameter). This value is indicative of the nebuliser's ability to produce particles suitable for a specific therapeutic treatment.

## 7 - DISPOSING OF THE APPARATUS AND ITS COMPONENTS

At the end of their useful life, the device and its components must be disposed of in accordance with the current regulations and environmental laws. If no legal obligations exist, the various materials must be sorted and disposed of separately.

The table below shows the materials of which the various components are made.

None of the materials used contains phthalates or traces of natural Latex.

COMPONENT	MATERIAL
Basic unit (1)	Waste from Electric and Electronic Equipment (WEEE)
External casing of the basic unit	ABS/Polystyrene
Nebuliser (9, 11)	Polycarbonate
Atomiser (10)	Nylon
Mask body (12) and fitting (13)	Silicone and polycarbonate
Mouthpiece (14)	Polycarbonate
Connecting tube (15)	Non-toxic PVC (phthalate-free polyvinyl chloride)
Box	Cardboard
Instructions	Paper

## 8 - ELECTROMAGNETIC COMPATIBILITY

This device complies with the requirements of IEC EN 60601-1-2:2015 on the electromagnetic compatibility of medical devices. Compliance with electromagnetic compatibility standards does not guarantee total immunity of the product. Some devices (mobile phones, pagers, etc.) may interrupt operation if used near medical devices.

### RECOMMENDED DISTANCES BETWEEN RADIO-COMMUNICATION DEVICES

The SOFFIO PLUS aerosol therapy device is designed to operate in an electromagnetic environment where radiated RF disturbance is under control. The customer or user of the SOFFIO PLUS aerosol therapy device can help prevent electromagnetic interference by leaving a minimum distance between mobile and portable RF communication devices (transmitters) and the SOFFIO PLUS aerosol therapy device. The table below shows the recommended distances for the SOFFIO PLUS aerosol therapy device based on the maximum output of the radio communication devices.

Nominal maximum transmitter output [W]	Separation distance [m] at the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \times \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \times \sqrt{P}$
0.01	0,12	0,12	0,23
0.1	0,38	0,38	0,73
1	1,2	1,2	2,33
10	3,8	3,8	7,3
100	12	12	23

The recommended separation distance  $d$  in metres [m] for transmitters with a nominal maximum transmitter output not shown in the table can be calculated using the equation applied to the transmitter frequency, where  $P$  is the nominal maximum transmitted output in Watt [W] declared by the manufacturer.

Note 1. For 80 MHz and 800 MHz, apply the distance for the higher frequency interval.

Note 2. These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and persons.

## 9 - OPERATING AND STORAGE CONDITIONS

### 9.1 - OPERATING CONDITIONS

- Temperature range +5 °C to +40 °C
- Relative humidity 15% to 93%
- Atmospheric pressure 700 hPa to 1060 hPa

### 9.2 - STORAGE CONDITIONS

- Temperature range -25 °C to +50 °C
- Relative humidity less than 93%
- Atmospheric pressure 700 hPa to 1060 hPa

## 10 - SPARE PARTS

Only original spare parts must be used.

## 11 - WARRANTY CONDITIONS

The manufacturer is liable for the safety, reliability and performance of the apparatus, provided that it is used in accordance with the instructions and for the intended purpose only, and any repairs are carried out by the manufacturer or an approved service centre. The manufacturer declines all liability and the warranty will be invalidated if the user fails to follow the instructions and the important notes contained herein.

## 12 - TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated accordingly to include any new or changed features.

If the booklet accompanying the apparatus gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

## 13 - COPYRIGHT

The information contained in this booklet refers only to the SOFFIO PLUS and to no other devices.

This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

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