
# Aidal



EN	Surgical aspirator	Operating instructions
IT	Aspiratore Chirurgico	Istruzioni d'uso
FR	Aspirateur chirurgical	Instructions d'utilisation
DE	Chirurgisches sauggerät	Bedienungsanleitung
ES	Aspirador Quirúrgico	Manual de instrucciones
PT	Aspirador cirúrgico	Instruções de uso
PL	Ssak Chirurgiczny	Instrukcje obsługi
NL	Chirurgische aspirator	Gebruiksaanwijzing
SV	Kirurgisk Aspirator	Bruksanvisning
FI	Kirurginen Aspirattori	Käyttöohjeet
EL	Αντλία Αναρροφήσης Για Χειρουργική Χρήση	Οδηγίες χρήσης
SK	Chirurgická Odsávačka	Návod na použitie
JA	医療用吸引器	使用説明書
KO	의료용 흡입기	사용안내서
CS	Kompaktní odsávačka	Návod k použití

# INDEX

EN	ENGLISH	p. 4
IT	ITALIANO	p. 7
FR	FRANÇAIS	p. 10
DE	DEUTSCH	p. 13
ES	ESPAÑOL	p. 16
PT	PORTUGUÊS	p. 20
PL	POLSKI	p. 23
NL	NEDERLANDS	p. 26
SV	SVENSKA	p. 29
FI	SUOMI	p. 32
EL	ΕΛΛΗΝΙΚΑ	p. 35
SK	SLOVENČINA	p. 38
JA	日本語	p. 41
KO	조선어	p. 44
CS	ČEŠTINA	p. 47

# Surgical Aspirator - Model: AIDAL and AIDAL DUAL BATTERY

Read these instructions carefully before using the device.

## 1 - INTENDED USE

The AIDAL surgical aspirator is a professional portable device intended for the suction of fluids and other substances from the human body. It can be used in the pharyngeal area and with tracheotomised patients. Vacuum pressure and flow can be adjusted. The AIDAL aspirator is not intended for thoracic drainage.

The product complies with the 'Restriction of use of certain hazardous substances in electrical and electronic equipment' of European Directive 2011/65/EU. On the basis of its intended use, the AIDAL aspirator is classified as a Class IIb medical device under the Medical Device Directive 93/42/EEC (and subsequent amendments) and is thus CE marked by the manufacturer, which is subject to surveillance by a Notified Body.

## 2 - ASPIRATOR DESCRIPTION

The aspirator is easy to handle and easy to use, robust, compact and quiet. It is equipped with a convenient carrying handle, comes with a reusable graduated collection jar and supports the use of disposable bags (sold separately) for up to 1000 ml of fluid. A safety stopper in the jar prevents aspirated fluids from penetrating into the aspirator unit if the maximum containable quantity is exceeded. The jar lid connects to an antibacterial filter that blocks the passage of solid particles, limits environmental contamination and blocks outflow of fluid in the event of a first fault.

The aspirator can be operated on a rechargeable 12V DC battery or a standard mains power supply.

AIDAL suction unit is also available with two 12VDC rechargeable batteries, identified on the label.

## 3 - PACKAGE CONTENTS

The package contains:

- The aspirator and rechargeable battery
- A power supply adapter and cord (cod.370021)
- A 1000 ml graduated collection jar
- Two antibacterial filters
- Tube kit (2 tubes of varying length)
- A reducing connector to connect a sterile disposable cannula\* and the connection tube
- Connector with valve for suction control
- One disposable cannula\*
- Operating instructions and warranty
- Aspirator carrying bag with strap

\* The cannula must comply with Directive 93/42 EEC concerning medical devices and must bear the CE marking.

\*\*Components 3 and 4 are not in contact with the patient

## 4 - OPTIONALS

- Cigarette lighter power supply socket 12V DC (only for recharge)
- 1000 ml disposable bags

NB: The operating instructions and intended use of the cigarette lighter socket are included in its relative packaging.

## 5 - OPERATING CONDITIONS

Maximum vacuum and maximum flow as indicated can be obtained in continuous operation under the following environmental conditions:

- Temperature 0 to 40 °C
- Relative humidity 20 to 93%
- Pressure 700 to 1060 hPa

## 6 - STORAGE CONDITIONS

The device can be stored under the following environmental conditions:

- Temperature -15 to +40 °C
- Relative humidity 0 to 93%
- Pressure 700 to 1060 hPa

⚠ Use appropriate personal protective equipment while operating the device.

## 7 - PREPARATION OF THE DEVICE (Figs. A, B, C)

### ► With collection jar

- a.1.1) Check the presence of the float (9,10,11) before installing the jar.
- a.1.2) Install the jar (1) in position.
- a.1.3) Fit the antibacterial filter (5), with the side marked "IN" towards the lid of the collection jar (2), into the lodging marked "VACUUM".
- a.1.4) Connect the short tube (6) to the aspirator unit socket (7) and antibacterial filter (5).
- a.1.5) Connect one end of the long tube (4) to the jar lid (2) connector fitting marked "PATIENT". Connect the other end of the tube to the cannula connector fitting (3) to which the disposable cannula is to be attached.
- a.1.6) Check that all components are correctly assembled.

### ► With disposable bag (sold separately)

- a.2.1) Insert the disposable bag into the collection jar (1).
- a.2.2) Connect one end of the first short tube (14) to the disposable bag (13) marked "VACUUM" and the other end to the antibacterial filter (5) with the marking "IN" towards the disposable bag.
- a.2.3) Connect one end of the second short tube (15) to the free side of the antibacterial filter and the other end of the tube to the aspirator unit socket (7).
- a.2.4) Connect one end of the long tube (16) to the disposable bag (13) connector fitting marked "PATIENT". Connect the other end of the tube to the cannula connector fitting (3) to which the disposable cannula is to be attached.
- a.2.5) Check that all components are correctly assembled.

⚠ The disposable bag is compatible with the collection jar supplied with the aspirator.

⚠ The vacuum system must always be used with the bacterial filter to prevent the contamination of the suction unit.

## 8 - DESCRIPTION OF CONTROLS & INDICATORS (Fig. D)

On the front of the aspirator are (Fig. D):

- b.1) ON/OFF green indicator light (17).
- b.2) Aspirator vacuum gauge with graduated kPa scale (18).
- b.3) Knob for adjustment of the vacuum minimum and maximum (19).
- b.4) Three indicator lights (green, yellow and red) to indicate the battery charge (20).
- b.5) Green indicator light to indicate the aspirator is connected to the mains power supply (21).

## 9 - OPERATING THE DEVICE ON MAINS POWER SUPPLY & BATTERY RECHARGING (Figs. D, E, F)

Fully charge the battery, before using the aspirator.

- c.1) Check that the mains power supply corresponds with that indicated on the power adapter.
- c.2) Connect the plug (24) of the power supply adapter to the socket of the aspirator (22), raising the water protection cover (23). Connect the adapter by the other end of the power supply cord (26) to the mains power.
- c.3) The green light (21) and an acoustic signal will indicate that the device is connected to the mains power.
- c.4) During charging the battery charge indicators (20) will indicate charge status. The green indicator will flash until the battery is fully charged (approximately 6 hours) and will then remain steadily lit.
- c.5) Press the power button (17) to turn on the aspirator.

⚠ Operating the device on mains power supply will also charge the battery.

⚠ Do not power on the aspirator if the power supply adapter or cord are damaged.

## 10 - OPERATING THE DEVICE ON INTERNAL BATTERY POWER SUPPLY (Fig. D)

Once the aspirator is in operation, the battery charge indicator lights (20) will indicate the battery charge status:

- The green indicator light indicates that the battery is fully charged.
- The yellow indicator light indicates that the battery is 30% charged.
- The red indicator light, together with a continuous acoustic signal, indicates that the battery is low and must be recharged as soon as possible.

Note = With the battery fully charged and the green indicator light (20) on the aspirator has about 25 minutes of charge under normal operation and with the yellow light on about 10 minutes. The red light indicates that the aspirator must immediately be connected to the mains power supply to ensure the continuity of patient care and recharging of the battery. Such operating times are indicative only and depend on the initial amount of battery charge. The acoustic signal continues to sound until the device is switched off or it is connected to the mains power supply.

⚠ Charge the battery at least once every 2 months when the device is not in use.

⚠ The battery has an average life span of 2 years. We recommend replacing it after that period.

## 11 - STARTING THE DEVICE & ADJUSTING THE VACUUM PRESSURE (Fig. D)

- d.1) Turn on the aspirator by pressing the green power button (17). Once activated, the button will light and the aspirator will start operation.

**Adjust the vacuum pressure by the vacuum adjustment dial (19) in the following way:**

- d.2) Close with a finger the open end of the cannula or tube.
- d.3) Turn the dial clockwise (19) to increase suction, anticlockwise to decrease it. Read the vacuum pressure value on the vacuum gauge.
- d.4) Once the desired value has been selected, remove your finger from the cannula or tube.
- d.5) Now the aspirator is ready for use.
- d.6) To turn off the device, press the green power button again (17).

N.B. = The kit contains a connector with valve that allows you to adjust the vacuum level without having to use the dial on the control panel of the device. To regulate the flow, open the valve and cover the hole gradually with your finger. The more you cover the hole, the higher the vacuum.

In order to avoid fluids entering the aspirator unit, empty the collection jar before the fluid reaches the float.

## 12 - DISPOSAL OF BODILY FLUIDS & DISPOSABLE BAG

It is essential to remember that fluids aspirated into the collection jar are potentially infectious and their disposal must be in accordance with procedures established for such waste.

Disposal of the jars, which are intended for the collection of potentially infectious body fluids, as well as disposal of the fluids themselves, must therefore be effectuated following hospital procedures for such types of waste. In the case of home use, prescriptions to use the aspirator must contain instructions for the disposal of single-use jars and aspirated liquids.

## 13 - CLEANING & DISINFECTION OF THE ASPIRATOR

The aspirator can be cleaned and disinfected using a sponge or soft cloth moistened with a cold disinfectant solution (e.g. a sodium hypochlorite solution with a concentration no greater than 2%), as directed by the manufacturer of the product. Ensure that the solution is compatible with the materials of the aspirator or follow hospital procedures.

- ⚠ Do not wash or immerse the aspirator under water.
- ⚠ The aspirator must always be turned off and disconnected from the mains while cleaning.

**14 - CLEANING & STERILIZATION OF REUSABLE COLLECTION JAR**

The jars must be cleaned and sterilized after each use, as described below:

**Cleaning**

Separate the main components of the jar such as the graduated vessel, the lid and the float. Wash all surfaces of the device thoroughly with hot water at a temperature not exceeding 60°C. Dry the parts thoroughly with a clean, soft non-abrasive cloth.

**Sterilization**

Place the parts in an autoclave for a sterilization cycle with a steam temperature of 121°C for at least 15 minutes (preferably not exceeding 30 minutes) taking care to position the graduated jar upside down (with its bottom upward).

Steam sterilization may be performed for a maximum of 15 cycles without adversely affecting the characteristics of the component materials used and therefore functionality. Sterilization may alter the colour or appearance of the components, but this will not affect operation of the device. However, a visual inspection of all the components is recommended before each assembly following sterilization. Sterilization can also be carried out cold by ethylene oxide or gamma rays following the instructions supplied with respective equipment.

**15 - ANTIBACTERIAL FILTER & CONNECTION TUBES**

The antibacterial filter and the PVC connection tubes are not steam-sterilizable and must be considered single-use. The bacterial filter must be replaced in case of contamination, yellowing or decrease in performance of the vacuum system. The filter must be replaced whenever the vacuum system is used on multiple patients. In any case, the filter should be replaced at least once every two months.

**16 - SAFETY FLOAT**

During cleaning and sterilization of the reusable jar, it is essential that the float inside the lid of the jar is removed and later reassembled for use. The operations to be performed are:

**Disassembly (Fig. B)**

- e.1) Open the jar, remove the float guide (9) from its lodging in the jar lid (12).
- e.2) Remove the float body (10) from the guide (9) and remove the rubber stopper.

**Assembly**

- e.3) Insert the rubber stopper (11) onto the float body (10) with the widest part of the stopper facing upward.
- e.4) Insert the float body (10) with the stopper (11) into the float guide (9).
- e.5) Insert the float guide (10) onto its lodging in the lid of the jar (12).
- e.6) Check that the float and stopper move freely inside the float guide.

**17 - DISPOSAL OF THE ASPIRATOR**

The aspirator and its components are made with plastic and metal materials. Their disposal must comply with all applicable regulations.

PART	MATERIAL
Base unit	Waste from electrical and electronic equipment (WEEE)
Battery	Lead
Outer shell of the base unit	ABS V0
Liquid collection tray	Polycarbonate
Tray lid	Polypropylene
Single-use tube	PVC
Cannula reducer	Nylon
Transport bag	Nylon
Box	Cardboard
Leaflet	Paper
Accessory bag and blue bag	Low density polyethylene (PE-LD)
Other accessories	See information on the packaging

**18 - PROBLEMS & POSSIBLE CAUSES & SOLUTIONS**

PROBLEM	POSSIBLE CAUSE	SOLUTION
• The device does not power on.	• Power button not activated. • Battery out of charge. • Power supply cord not connected.	• Press the power button.  • Recharge the battery. • Connect the power supply cord.
• No suction.	• Connector tubes or antibacterial filter not installed properly. • Jar lid not closed properly. • Disposable bag not fitted properly in the collection jar. • Vacuum adjustment dial rotated fully anticlockwise. • Jar float stuck.	• Install tubes and filter as indicated in the instructions.  • Close the jar lid properly.  • Fit the disposable bag properly.  • Rotate the dial clockwise.  • Free the float and check that it can move freely in its lodging.

PROBLEM	POSSIBLE CAUSE	SOLUTION
• The performance of the device is different from that stated.	• Low battery charge. • Power adapter not properly connected to mains or aspirator. • Vacuum adjustment dial rotated fully anticlockwise. • Jar lid not closed properly. • Disposable bag not fitted properly in the collection jar. • Jar float stuck.  • Connection tubes faulty.	• Recharge the battery. • Connect the power adapter properly.  • Rotate the dial clockwise.  • Close the jar lid properly.  • Fit the disposable bag properly.  • Free the float and check that it can move freely in its lodging. • Replace the tubes.

If, after checking the above, the aspirator does not work properly, inspection at an authorized laboratory is necessary.

The following maintenance operations must be carried out by qualified technical assistance staff at an approved laboratory, at the intervals specified:

- replacement of the battery whenever the life of the latter is not optimal or at least every two years.
- replacement of air filter, each time the performance of the suction unit (vacuum rate and flow rate) is lower than the rated value or at least every two years.

**19 - TECHNICAL DATA**

	AIDAL	AIDAL DUAL BATTERY (version with two batteries)
Power supply adapter	output 15V DC 4.0A 60W	
Voltage AIDAL	12V DC	
Definition according to UNI EN ISO 10079-1	high vacuum / High flow	
Consumption at 12V DC	36 W	
Maximum vacuum	-87 kPa (-0.87 bar) ±8%	
Maximum air flow	30±4 l/min	
Battery charging time	7 hours approx.	
Battery type	Pb 12V DC	2 x Pb 12V DC
Battery powered operating time	40 minutes (depending on use)	60 minutes (depending on use)
Type of operation	Continuous	
Noise level pursuant to UNI EN 10079-1:2009	55,5 dBA	
Standard jar capacity	1000ml	
Weight with 1000ml jar	3,8 kg	4,6 kg
Dimensions: length width height	350mm 140mm 280mm	
Classification according to CEI EN 60601-1	<ul style="list-style-type: none"> <li>• Class II device.</li> <li>• Device with internal power source.</li> <li>• Device with Type BF applied part.</li> <li>• Device not suitable for use with anaesthetic mixtures with air, oxygen or nitrous oxide.</li> </ul>	
Classification according to UNI EN ISO 10079-1	<ul style="list-style-type: none"> <li>• Device for high vacuum and high flow.</li> <li>• Device suitable for use in the field and in ambulances.</li> </ul>	

**20 - ELECTROMAGNETIC COMPATIBILITY**

This device complies with the requirements of CEI EN 60601-1-2:2015 for the electromagnetic compatibility of medical devices. Compliance with EMC standards does not guarantee the total immunity of the product; some devices (cell phones, pagers, etc.), when used near medical equipment, can stop it from functioning properly.

ELECTROMAGNETIC EMISSION		
Aidal is suitable for use in the specified electromagnetic environment. The customer or the user of the Aidal should assure that it is used in an electromagnetic environment as described below:		
Emissions test	Compliance	Electromagnetic Environment
RF emissions CISPR 11	Group 1	Aidal aspirator uses RF energy only for its internal function. Therefore, its RF emission is very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Aidal is suitable for use in all establishments, including domestic establishment and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions EN 61000-3-2	Class A conforming	
Voltage fluctuations/ flicker emissions EN 61000-3-3	Complies	

ELECTROMAGNETIC IMMUNITY			
Aidal is intended for use in the electromagnetic environment specified below. The customer or user should assure that it is used in an electromagnetic environment.			
Immunity Test	EN 60601-1-2 Test level	Compliance level	Electromagnetic Environment - guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst EN 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Immunity Test	EN 60601-1-2 Test level	Compliance level	Electromagnetic Environment - guidance
Surge EN 61000-4-5	± 2 kV differential mode ± 1 kV common mode	± 2 kV differential mode ± 1 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines.  EN 61000-4-11	<5 % UT for 0,5 cycle (>95 % dip in UT)	<5 % UT for 0,5 cycle (>95 % dip in UT)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Aidal requires continued operation during power mains interruptions, it is recommended that the Aidal be powered from an uninterruptible power supply or a battery.
	40 % UT for 5 cycles (60 % dip in UT)	40 % UT for 5 cycles (60 % dip in UT)	
	70 % UT for 25 cycles (30 % dip in UT)	70 % UT for 25 cycles (30 % dip in UT)	
	<5 % UT for 5 s (>95 % dip in UT)	<5 % UT for 5 s (>95 % dip in UT)	
Power frequency (50/60 Hz) magnetic field EN 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity Test	EN 60601-1-2 Test level	Compliance level	Electromagnetic Environment
Radiated RF EN 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of Aidal aspirator, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter  <b>Recommended separation distance</b> $d = 1.2 \times \sqrt{P}$ 80 MHz to 800MHz $d = 2.3 \times \sqrt{P}$ 800 MHz to 2.5GHz $d = 1.2 \times \sqrt{P}$ Where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths for fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF EN 61000-4-6	3 V 150 kHz to 80 MHz	3 V	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structure, objects and people.			

RECOMMENDED DISTANCES FROM RADIO-COMMUNICATION DEVICES			
The AIDAL surgical aspirator is intended for use in an electromagnetic environment in which radiated RF disturbances are under control. The customer or the operator of the AIDAL device can help prevent electromagnetic interference by maintaining a minimum distance between mobile and portable RF communication equipment (transmitters) and the AIDAL device as recommended below in relation to the peak output of the radio-communication devices.			
Maximum rated power output of the transmitter [W]	Separation distance at the frequency of the transmitter [m]		
	150 kHz - 80 MHz $d = 1.2 \times \sqrt{P}$	80 MHz - 800 MHz $d = 1.2 \times \sqrt{P}$	800 MHz - 2.5 GHz $d = 2.3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters with a maximum rated output power not listed above, the recommended separation distance d in metres (m) can be calculated using the equation applicable to the frequency of the transmitter, where P is the maximum rated output power of the transmitter in Watts (W) according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz the separation distance is applied for the higher frequency range. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

**21 - WARNINGS**

The product has been designed with particular attention to patient and user safety. Nevertheless, the following precautions must be taken.

- The aspirator must only be used by trained and qualified personnel.
- Do not modify this equipment without the manufacturer's authorization.
- Before operating the aspirator in any way, read all the instructions carefully, paying particular attention to warnings and how to prepare and use the device.
- Check that every precaution has been taken to avoid any danger arising from accidental contact with blood or bodily fluids.
- The aspirator must only be used for suction as indicated in this operating instruction manual as in the paragraph "Intended Use". Other uses are not permitted.
- The aspirator is not suitable for thoracic drainage.
- The vacuum generated in the aspirator suction outlet, in its connecting tubes or collection jar may represent a potential hazard.

- The connecting tubes (PVC) between the various components are for single-use.
- If the antibacterial filter and the safety float do not intervene and aspirated fluid gets into the aspirator unit, immediately interrupt the power supply.
- The device is protected against access to hazardous parts and against vertically falling drops of water (IP34). IP protection is no longer guaranteed while the unit is being recharged.
- The aspirator is not suitable for use in the presence of anaesthetic mixtures with air, oxygen or nitrous oxide.
- Do not touch live parts such as the power cord, power supply adapter and the power button with wet hands.
- In case of an accidental fall into water the aspirator unit can be extracted only after interrupting the power supply. After such an occurrence the suction unit may not be used and requires a complete overhaul.
- In use, the aspirator must be kept on a hard level surface to prevent leakage of fluids from the collection jar. The device must not be used in movement.
- Ensure that the lid of the jar is properly closed in order to avoid accidental spillage of fluids.
- In use, the aspirator must be kept away from source of heat.
- The power cord and connecting tubes, for their length, may pose a risk of strangulation.
- Some components for their small size, if swallowed, may present a risk of suffocation.
- It is recommended to fully charge the aspirator before each use.
- In case of allergic reaction due to contact with the materials of this device, seek medical attention.
- The manufacturer has defined a useful lifespan of at least five years from first use.
- The battery has an average life span of 2 years. We recommend replacing it after that period.

## 22 - WARRANTY

See the warranty certificate inside the package.

Accessories and batteries are covered by warranty against workmanship defects. The manufacturer is liable for the safety, the reliability and the performance of the device if the same is used in accordance with the operating instructions and its intended use. The manufacturer is liable for the safety, the reliability and the performance of the device if repairs are carried out at the manufacturer or at a laboratory authorised by the same. If the user fails to comply with the Operating Instructions and Warnings contained in this manual, the manufacturer's liability and warranty conditions shall be considered void.

## 23 - COPYRIGHT

All information contained herein may not be used for purposes other than originally intended. This manual is owned by Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part without the prior written permission of Air Liquide Medical Systems S.r.l. All rights reserved.

## 24 - UPDATES TO TECHNICAL CHARACTERISTICS

In order to continuously improve performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are periodically reviewed and improved. Operating instruction manuals are therefore amended to ensure constant compliance with the characteristics of the devices launched to market. If the instruction manual accompanying this device is lost, a copy corresponding to the supplied device may be obtained from the manufacturer by specifying the data shown on the product identification label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.