Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Brugsanvisning
Instruksjoner for bruk / Användningsinstruktioner / Käyttöohjeet / Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Inštrukcie na používanie / Használati útmutató
Instruțiuni de folosire / Navodila za uporabo / Οδηγίες χρήσης / Инструкции по применению / Інструкція із застосування / 使用説明 / 使用説明書 / 사용 지침


Alnest N1

Masque nasal avec fuites d'expiration / Nasal mask with exhalation orifices / Nasenmaske mit Ausatemungsöffnungen / Mascarilla nasal con orificios de exhalación / Maschera nasale con fori di esalazione / Máscara nasal com orificios de exalação / Nasalmaske med udåndingsåbninger / Nesemaske som er uten utpustningsåpninger / Näsmask med utandningshål / Nenänaamari jossa on uloshengityksaukot / Neusmasker met uitademingsopeningen / Maska nosowa z otworami wentylacyjnymi / Nosní maska s výdechovými otvory / Nazálna maska s výdychovými otvormi / Orrmaskz kilégzőnyílásokkal / Mască nazală cu orificii de dispersare a aerului expirat / Nosna maska s luknjami za izdihavanje / Ρινική μάσκα με οπές εκπνοής / Назальная маска с вентиляционными отверстиями / Назальная маска з вентиляційними отворами / 带有标准通气口的鼻面罩 / 鼻マスク通気穴付き / 호기구멍이 있는 비강 마스크서


A


B


①

②


③


B


④











⑤

⑥

C**1****2****D****1****3****2****4**

E**1****2****3****4**

FRANÇAIS	p. 09	ROMÂNĂ.....	p. 58
ENGLISH.....	p. 12	SLOVENŠČINA.....	p. 61
DEUTSCH	p. 15	ΕΛΛΗΝΙΚΑ	p. 64
ESPAÑOL	p. 19	РУССКИЙ.....	p. 67
ITALIANO	p. 22	УКРАЇНЕЦЬ.....	p. 70
PORTUGUÊS	p. 25	繁體中文.....	p. 74
DANSK.....	p. 28	日本語.....	p. 76
NORSK	p. 31	한국어.....	p. 79
SVENSKA.....	p. 35		
SUOMI.....	p. 38		
NEDERLANDS.....	p. 41		
POLSKI	p. 45		
ČEŠTINA.....	p. 48		
SLOVENČINA.....	p. 51		
MAGYAR.....	p. 54		

1 - UTILISATION

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1 - INTENDED USE

ALNEST N1 is a reusable nasal mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by multiple adult patients (weighing more than 30 kg) who have been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

WARNINGS

- The mask should be used exclusively with positive pressure equipment recommended by your doctor or respiratory therapist supplying a minimum pressure of 4 cmH₂O: at lower pressures, exhaled air safety evacuation through the exhalation orifices is not guaranteed and so a partial re-inhalation may occur.
- Before using the mask, it needs to check that the exhalation orifices are not occluded, even partially, to not jeopardize the safety and the quality of the therapy.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual, or exceeding the recommended number of regeneration cycles, can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2 - MASK COMPONENTS (section A)

a) cushion; b) frame; c) elbow connector; d) swivel fitting; e) headgear; f) buckle.

3 - ASSEMBLING AND WEARING THE MASK (section B)

- Bring the gray cushion (a) ring closer to the frame (b) so as to first couple the upper part of the cushion to the upper tooth of the frame, then place your fingers and push to also hook the lower part of the cushion to the frame (B.1).
- Holding the mask gently in close contact to the face, wear the headgear (e) over the head making sure the cushion (a) adheres to the profile around the nose (B.2).
- Attach the two lower straps of the headgear (e) to the hooks of the frame (b) so that they pass under the ears (B.3).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.4.
- Connect the swivel fitting (d) to the equipment tube (B.5), then turn on the equipment at the prescribed pressure.
- Lie down and, if necessary, adjust the headgear again (B.6).

4 - USING THE MASK

To properly remove, disassemble and reassemble the mask, refer to the respective illustrations in section C - D - E.

5 - CLEANING THE MASK

For daily cleaning, the mask must be disassembled in its components, each of which must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

To clean the headgear, perform the same procedure weekly.

For a deeper decontamination of the mask, consult the “Guide to decontamination, disinfection and sterilization” available electronically on the website medicaldevice.airliquide.com or at your dealer, together with the instructions provided by the manufacturer of the disinfectant solution, especially regarding the compatibility with the materials of the various components.

⚠ PRECAUTIONS

- Remove the mask after disconnecting at least one of the lower headgear straps.
- Do not open the straps to avoid having to adjust the headgear again the next time you use it.
- Be very careful to remove any dirt from the exhalation holes.
- The headgear can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

6 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS

The mask can be regenerated for up to 20 cycles by one of the following processes:

- high level thermal disinfection;
- high level chemical disinfection.

For full detailed mask regeneration instructions, consult the “Guide to decontamination, disinfection and sterilization” available electronically on the website medicaldevice.airliquide.com or at your dealer.

⚠ PRECAUTIONS

- *The mask must always be regenerated before use on a new patient.*
- *The headgear must be replaced.*

7 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	<ul style="list-style-type: none"> - Flow generator turned off or not working. - Breathing circuit ventilator not inserted correctly. 	<ul style="list-style-type: none"> - Switch on or adjust the flow generator. - Connect the breathing circuit ventilator to the flow generator and the mask.
Noisy and annoying air leaks.	<ul style="list-style-type: none"> - The mask has not been assembled correctly. - The mask is not positioned correctly. - The size of the mask may not be correct. 	<ul style="list-style-type: none"> - Disassemble the mask and then reassemble it as described in the user manual. - Reposition the mask on your face and adjust the headgear straps as described in the user manual. - Contact your doctor or therapy technician to find the right size.
The exhalation orifices become dirty to the point of occlusion.	Incorrect or insufficient daily cleaning.	Clean the exhalation orifices using a soft-bristled toothbrush.


8 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + PBT
Frame (b)	Polycarbonate
Elbow connector (c)	Polycarbonate + Nylon
Swivel fitting (d)	Copolyester
Headgear (e)	Nylon/UBL/Polyurethane
Buckle (f)	Nylon

9 - TECHNICAL DATA

PRESSURE - FLOW CURVE*

*Due to machining differences, the outflow from the exhalation orifices may vary ($\pm 10\%$).

Regulatory information	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 17510
Operating pressure range	4 cmH ₂ O / 32 cmH ₂ O
Operating temperature range	+5 °C / +40 °C
Dead space	Size S: 64 ml Size M: 80 ml Size L: 96 ml
Resistance	at 50 l/min: 0,24 cmH ₂ O at 100 l/min: 0,93 cmH ₂ O
A-Weighted Sound Pressure (UNI EN ISO 3744)	16,7 ± 3,2 dB(A)
A-Weighted Sound Power (UNI EN ISO 3744)	24,7 ± 3,2 dB(A)

Dimensions (without headgear)	Size S: 150 x 80 x 90 mm Size M: 150 x 80 x 90 mm Size L: 150 x 80 x 95 mm
Weight	Size S: 83 g Size M: 84 g Size L: 85 g

10 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components. The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market.

If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.