Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Brugsanvisning / Instruksjoner for bruk / Användningsinstruktioner / Käyttöohjeet / Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Inštrukcie na používanie / Használati útmutató / Instruțiuni de folosire / Navodila za uporabo / Οδηγίες χρήσης / Инструкции по применению / Інструкція із застосування / 使用說明 / 使用說明書 / 사용 지침 / تعليمات الاستخدام / הוראות שימוש




# Alnest **SOFT**



## Baby





### Non Vented

Masque nasal / Nasal mask / Nasenmaske / Mascarilla nasal / Maschera nasale / Máscara nasal / Nasalmaske / Nesemaske / Näsmask / Nenänaamari / Neusmasker / Maska nosowa / Nosní maska / Nazálna maska / Orrmaszk / Mască nazală / Nosna maska / Ρινική μάσκα / Назальная маска / Назальна маска / 鼻面罩 / 鼻マスク / 비강 마스크서 / قناع أنفي / סוכת נא




**OPTIONAL**


1


2



3

4

FRANÇAIS .....	p. 08	SLOVENČINA.....	p. 49
ENGLISH.....	p. 11	MAGYAR.....	p. 52
DEUTSCH .....	p. 14	ROMÂNĂ.....	p. 55
ESPAÑOL .....	p. 17	SLOVENŠČINA .....	p. 58
ITALIANO .....	p. 21	ΕΛΛΗΝΙΚΑ .....	p. 61
PORTUGUÊS .....	p. 24	РУССКИЙ.....	p. 64
DANSK.....	p. 27	УКРАЇНСЬКА.....	p. 68
NORSK .....	p. 30	繁體中文 .....	p. 71
SVENSKA.....	p. 33	日本語 .....	p. 73
SUOMI.....	p. 36	한국어.....	p. 76
NEDERLANDS.....	p. 39	العربية.....	p. 81
POLSKI .....	p. 43	עברית.....	p. 85
ČEŠTINA .....	p. 46		

## 1 - USAGE PRÉVU

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## ENGLISH

## 1 - INTENDED USE

This medical device belongs to the category "Masks for non-invasive ventilation and related connectors", that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

## INDICATION FOR USE AND INTENDED USERS

Alnest SOFT Baby Non Vented is a reusable nasal mask without calibrated exhalation orifices intended to be used at home or in a hospital setting by multiple pediatric patients (weighing between 3 and 12 kg) in controlled non-invasive positive pressure ventilation (NIPPV).

 WARNINGS

- The mask must be used under the supervision of an adult who has read and understands the instructions.
- The mask should only be used with a ventilator fitted with an active exhalation valve and/or with a respiratory circuit with an exhalation valve recommended by your doctor or respiratory therapist.

- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual, or exceeding the recommended number of regeneration cycles, can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- To ensure proper operation of the hose assembly, handle it with care, especially avoiding pulling it.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## 2 - MASK COMPONENTS (section A)

- a) cushion; b) hose assembly; c) sampling port cap; d) headgear.

**3 - ASSEMBLING AND WEARING THE MASK (section B)**

- Connect the hose assembly (b) by inserting the sphere on the elbow fitting in the rigid slot in the cushion (a) (B.1).
- Insert the three top straps of the headgear (d) into the three top slots of the pad (a), then fasten the Velcro straps (B.2).
- Holding the mask gently in close contact to the face, wear the headgear (d) over the child's head making sure the cushion (a) adheres to the profile around the nose (B.3).
- Insert the two bottom straps into the two bottom slots of the pad (a), then fasten the Velcro straps (B.4).
- Slightly pull the straps to ensure the stability of the mask according to the sequence B.5.
- Connect the blue fitting of the hose assembly (b) to the equipment tube (B.6), then turn on the equipment at the prescribed pressure.

NOTE: Alternatively, the mask assembly can also be completed by removing the forehead rest from the cushion as illustrated in section OPTIONAL C.

**4 - DISASSEMBLING THE MASK (section D)**

To properly disassemble the mask, refer to the illustrations in section D.

**5 - CLEANING THE MASK**

For daily cleaning, the mask must be disassembled in its components, each of which (except for the headgear) must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

**To clean the headgear, perform the same procedure weekly.**

For a deeper decontamination of the mask, consult the "Guide to decontamination, disinfection and sterilization" available electronically on the website [it.medicaldevice.airliquide.com](http://it.medicaldevice.airliquide.com) or at your dealer, together with the instructions provided by the manufacturer of the disinfectant solution, especially regarding the compatibility with the materials of the various components.

**⚠ PRECAUTIONS**

- Remove the mask after disconnecting at least one of the lower headgear straps.
- The headgear can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

**6 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS**

The mask can be regenerated for up to 12 cycles by one of the following processes:

- high level thermal disinfection;
- high level chemical disinfection;
- autoclave sterilization at 121°C for a maximum of 30 minutes.

For full detailed mask regeneration instructions, consult the "Guide to decontamination, disinfection and sterilization" available electronically on the website [it.medicaldevice.airliquide.com](http://it.medicaldevice.airliquide.com) or at your dealer.

**⚠ PRECAUTIONS**

- *The mask must always be regenerated before use on a new patient.*
- *The headgear and hose assembly must be replaced.*

## 7 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	<ul style="list-style-type: none"> <li>- Flow generator turned off or not working.</li> <li>- Breathing circuit ventilator not inserted correctly.</li> </ul>	<ul style="list-style-type: none"> <li>- Switch on or adjust the flow generator.</li> <li>- Connect the breathing circuit ventilator to the flow generator and the mask.</li> </ul>
Noisy and annoying air leaks.	<ul style="list-style-type: none"> <li>- The mask has not been assembled correctly.</li> <li>- The mask is not positioned correctly.</li> <li>- The size of the mask may not be correct.</li> </ul>	<ul style="list-style-type: none"> <li>- Disassemble the mask and then reassemble it as described in the user manual.</li> <li>- Reposition the mask on your face and adjust the headgear straps as described in the user manual.</li> <li>- Contact your doctor or therapy technician to find the right size.</li> </ul>

## 8 - DISPOSAL

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + Polycarbonate
Hose assembly (b)	Polycarbonate + Ethylene-Octene Copolymer + POM
Sampling port cap (c)	Silicone
Headgear (d)	Nylon/UBL/Polyurethane

## 9 - TECHNICAL DATA

<b>Regulatory information</b>	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 5356-1
<b>Operating pressure range</b>	4 cmH <sub>2</sub> O / 24 cmH <sub>2</sub> O
<b>Operating temperature range</b>	+5 °C / +40 °C
<b>Dead space</b>	Size XXS: 39.8 ml Size XS: 39.8 ml Size S: 39.8 ml
<b>Dimensions (without headgear)</b>	Size XXS: 90 x 90 x 50 mm Size XS: 90 x 90 x 50 mm Size S: 90 x 90 x 50 mm
<b>Weight</b>	Size XXS: 48 g Size XS: 50 g Size S: 52 g

## 10 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer. For more information on the rights under the warranty, please contact your retailer.

**COPYRIGHT**

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

**TECHNICAL CHARACTERISTICS UPDATES**

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market.

If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.