
Mode d'emploi / Instructions for use / Gebrauchsanweisung / Instrucciones de uso / Istruzioni per l'uso / Instruções de uso / Brugsanvisning / Instruksjoner for bruk / Användningsinstruktioner / Käyttöohjeet / Gebruiksaanwijzing / Instrukcja użycia / Návod k použití / Návod na použitie / Használati útmutató / Instrucțiuni de utilizare / Navodila za uporabo / Οδηγίες χρήσης / Инструкция по использованию / Інструкція із застосування / 使用说明 / 使用說明書 / 사용 지침 / تعليمات الاستخدام / תורות שימוש







# Alnest **SOFT**





## Baby

Masque nasal avec fuites d'expiration / Nasal mask with exhalation orifices / Nasenmaske mit Ausatemöffnungen / Mascarilla nasal con orificios de exhalación / Maschera nasale con fori di esalazione / Máscara nasal com orifícios de exalação / Nasalmaske med udåndingsåbninger / Nesemaske med utpustningsåpninger / Näsmask med utandningshål / Nenämaski uloshengitysaukoilla / Neusmasker met uitademingsopeningen / Maska nosowa z otworami wydechowymi / Nosní maska s výdechovými otvory / Nazálna maska z výdychovými otvormi / Orrmaszk kilégzőnyílásokkal / Mască nazală cu orificii de expirație / Nosna maska z luknjami za izdihavanje / Ρινική μάσκα με οπές εκπνοής / Назальная маска с вентиляционными отверстиями / Назальна маска з вентиляційними отворами / 带排气口的鼻罩 / 呼吸排出孔付き鼻マスク / 호기구가 있는 비강 마스크 / قناع أنفي مع فتحات للزفير / قناع أنفي עם חרירי נשיפה / מסכת אף-פה עם חרירי נשיפה



OPTIONAL


1


2




3

4

## INDEX

FRANÇAIS .....	8	ČEŠTINA .....	51
ENGLISH .....	11	SLOVENČINA .....	54
DEUTSCH .....	15	MAGYAR .....	58
ESPAÑOL .....	18	ROMÂNĂ .....	61
ITALIANO .....	22	SLOVENŠČINA .....	65
PORTUGUÊS .....	26	ΕΛΛΗΝΙΚΑ .....	68
DANSK .....	29	РУССКИЙ .....	72
NORSK .....	33	УКРАЇНСЬКА .....	76
SVENSKA .....	36	简体中文 .....	80
SUOMI .....	40	日本語 .....	83
NEDERLANDS .....	43	한국어 .....	86
POLSKI .....	47	العربية .....	93
		עברית .....	97

## 1 - USAGE PRÉVU

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## 1 - INTENDED USE

This medical device belongs to the category "Masks for non-invasive ventilation and related connectors", that includes nasal, nasal pillow and oro-nasal masks, along with the related connectors, which are intended to be used at home or in a hospital setting by a single or by multiple patients who need for a non-invasive positive pressure ventilation (NIPPV) therapy.

## INDICATION FOR USE AND INTENDED USERS

Alnest SOFT Baby is a reusable nasal mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by multiple pediatric patients (weighing between 3 and 12 kg) who have been prescribed for a non-invasive positive pressure ventilation therapy (NIPPV), e.g. continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

**⚠ WARNINGS**

- The mask must be used under the supervision of an adult who has read and understands the instructions.
- The mask should be used exclusively with positive pressure equipment recommended by your doctor or respiratory therapist supplying a minimum pressure of 4 cmH<sub>2</sub>O: at lower pressures, exhaled air safety evacuation through the exhalation orifices is not guaranteed and so a partial re-inhalation may occur.
- Before using the mask, it needs to check that the exhalation orifices are not occluded, even partially, to not jeopardize the safety and the quality of the therapy.
- In case of discomfort, irritation or allergic reactions to any component of the mask, consult your doctor or respiratory physiotherapist.
- In case of additional oxygen administration, smoking or using open flames is prohibited.
- When using oxygen and the equipment is not in operation, turn off the oxygen dispenser.
- Do not use the mask in case of vomiting or nausea.
- Keep away from light.
- It is recommended to wash the mask before the first use.
- The mask shall be cleaned after each use.
- Frequency of cleaning, methods of cleaning or the use of cleaning agents other than those specified in this user manual, or exceeding the recommended number of regeneration cycles, can jeopardize the safety and the quality of the therapy.
- Before use, check the integrity of the mask. If damaged during transport, notify your dealer.
- In case of visible deterioration (cracks, tears, etc.), it needs to eliminate and replace the damaged mask component.
- Do not leave the components of the mask unattended, some of these may be swallowed by children.
- To ensure proper operation of the hose assembly, handle it with care, especially avoiding pulling it.
- Contact your dealer as a healthcare professional in case of functional problems, size problems or questions regarding the application of the medical device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**2 - MASK COMPONENTS (section A)**

a) pad; b) hose assembly; c) swivel fitting; d) sampling port cap; e) headgear.

**3 - ASSEMBLING AND WEARING THE MASK (section B)**

- Connect the hose assembly (b) by inserting the sphere on the elbow fitting in the rigid slot in the cushion (a) (B.1).
- Connect the large part of the swivel fitting (c) to the free end of the hose assembly (b) (B.2).
- Insert the three top straps of the headgear (e) into the three top slots of the pad (a), then fasten the Velcro straps (B.3).
- Holding the mask gently against the face, slip the headgear (e) over the head of the child, making sure that the pad (a) adheres to the profile around the nose (B.4).
- Insert the two bottom straps into the two bottom slots of the pad (a), then fasten the Velcro straps (B.5).
- Slightly pull the straps to ensure the mask is firmly in place according to the sequence in B.6.
- Connect the swivel fitting (c) of the hose assembly to the device hose (B.7), then turn on the device at the prescribed pressure.

NOTE: Alternatively, the mask assembly can also be completed by removing the forehead rest from the cushion as illustrated in section OPTIONAL C.

**4 - DISASSEMBLING THE MASK (section D)**

To properly disassemble the mask, refer to the illustrations in section D.

**5 - CLEANING THE MASK**

For daily cleaning, the mask must be disassembled in its components, each of which (except for the headgear) must be carefully washed in warm water (about 30°C) using neutral soap and then rinsed thoroughly with clean fresh water. Finally, the parts must be air dried, away from direct sunlight.

**To clean the headgear, perform the same procedure weekly.**

For a deeper decontamination of the mask, consult the "Guide to decontamination, disinfection and sterilization" available electronically on the website [it.medicaldevice.airliquide.com](http://it.medicaldevice.airliquide.com) or at your dealer, together with the instructions provided by the manufacturer of the disinfectant solution, especially regarding the compatibility with the materials of the various components.

**⚠ PRECAUTIONS**

- Remove the mask after disconnecting at least one of the lower headgear straps.
- Be very careful to remove any dirt from the exhalation holes.
- The headgear can only be washed in water.
- Do not clean the mask with solutions containing alcohol, aromatic compounds, humidifiers, antibacterial agents, hydrating agents or aromatic oils as they would damage the product and reduce its life.
- Do not leave the mask immersed for more than 10 minutes.
- Do not clean the mask in the dishwasher.
- Check carefully all the components to make sure there are no soap deposits or foreign bodies inside the mask.

**6 - REGENERATION OF THE MASK FOR REUSE ON DIFFERENT PATIENTS**

The mask can be regenerated for up to 12 cycles by one of the following processes:

- high level thermal disinfection;
- high level chemical disinfection;
- autoclave sterilization at 121°C for a maximum of 30 minutes.

For full detailed mask regeneration instructions, consult the "Guide to decontamination, disinfection and sterilization" available electronically on the website [it.medicaldevice.airliquide.com](http://it.medicaldevice.airliquide.com) or at your dealer.

**⚠ PRECAUTIONS**

- *The mask must always be regenerated before use on a new patient.*
- *The headgear and hose assembly must be replaced.*

**7 - TROUBLESHOOTING**

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The airflow doesn't reach the mask.	- Flow generator turned off or not working. - Breathing circuit ventilator not inserted correctly.	- Switch on or adjust the flow generator. - Connect the breathing circuit ventilator to the flow generator and the mask.
Noisy and annoying air leaks.	- The mask has not been assembled correctly. - The mask is not positioned correctly. - The size of the mask may not be correct.	- Disassemble the mask and then reassemble it as described in the user manual. - Reposition the mask on your face and adjust the headgear straps as described in the user manual. - Contact your doctor or therapy technician to find the right size.
The exhalation orifices become dirty to the point of occlusion.	Incorrect or insufficient daily cleaning.	Clean the exhalation orifices using a soft-bristled toothbrush.


**8 - DISPOSAL**

The mask components are made of the materials indicated below. They do not contain hazardous substances and can therefore be disposed of with ordinary household waste.

COMPONENT	MATERIAL
Cushion (a)	Silicone + Polycarbonate
Hose assembly (b)	Polycarbonate + TPE + Polypropylene
Swivel fitting (c)	Copolyester
Sampling port cap (d)	Silicone
Headgear (e)	Nylon/UBL/Polyurethane

## 9 - TECHNICAL DATA

### PRESSURE-FLOW CURVE

*\*Due to machining differences, the outflow from the exhalation orifices may vary ( $\pm 10\%$ ).*

<b>Regulatory information</b>	CE marking in accordance with Regulation 2017/745 (EU) and subsequent implementation amendments. Class IIa Main product standard UNI EN ISO 17510
<b>Operating pressure range</b>	4 cmH <sub>2</sub> O / 24 cmH <sub>2</sub> O
<b>Operating temperature range</b>	+5°C / +40°C
<b>Dead space</b>	Size XXS: 53.4 ml Size XS: 53.4 ml Size S: 53.4 ml
<b>Resistance</b>	at 50 l/min: 0.55 cmH <sub>2</sub> O at 100 l/min: 1.80 cmH <sub>2</sub> O

<b>A-Weighted Sound Pressure (UNI EN ISO 3744)</b>	15.8 $\pm$ 3.2 dB(A)
<b>A-Weighted Sound Power (UNI EN ISO 3744)</b>	19.5 $\pm$ 3.2 dB(A)
<b>Dimensions (without headgear)</b>	Size XXS: 90 x 90 x 50 mm Size XS: 90 x 90 x 50 mm Size S: 90 x 90 x 50 mm
<b>Weight</b>	Size XXS: 49 g Size XS: 51 g Size S: 53 g

## 10 - WARRANTY CERTIFICATE

Air Liquide Medical Systems S.r.l. guarantees its product against material or manufacturing defects for 180 days from the date of purchase, provided that the conditions of use in the instructions for use are observed.

If the product is faulty under normal use, Air Liquide Medical Systems S.r.l. will, at its discretion, repair or replace the faulty product or its components.

The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the operating instructions and its intended use.

Keep the original casing for the entire warranty period. The warranty in case of product defects can be applied by contacting your dealer.

For more information on the rights under the warranty, please contact your retailer.

### COPYRIGHT

All the information in this manual may not be used for purposes other than the original. This manual is the property of Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without written permission from the company. All rights are reserved.

### TECHNICAL CHARACTERISTICS UPDATES

For a continuous improvement of performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are regularly subjected to review and change.

Instruction manuals are therefore modified to ensure their constant consistency with the characteristics of devices placed on the market. If the instruction manual accompanying this device is lost, contact the manufacturer to get a copy of the version corresponding to the device by quoting the references shown on the label.

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AIR LIQUIDE MEDICAL SYSTEMS S.r.l.