

Respireo 

Respireo
Primo 



Easy to use, and comfortable, Respireo™ N and Respireo™ Primo N are adjustable **nasal masks** suitable for patients on CPAP and NIV.


Made in Italy, these masks provide quality design and comfort for your patients on a daily basis.

Simple & comfortable

 **Air Liquide**
HEALTHCARE

Respireo^N

Respireo^N Primo^N



• The silicone **forehead pad** provides stability while relieving the pressure point on the nasal bridge



• An **adjustment knob** designed to adjust the support on the patient's forehead (16 possible positions)

• The silicone **double-wall cushion** provides a comfortable and optimal seal by adapting to the patient's face



• Both **masks** are made of only 3 transparent parts (the frame, the elbow and the cushion), which makes it easy to assemble and disassemble for cleaning



• The **hooks** system makes the mask easy to put on and quick to remove in case of emergency

A 4-point connection headgear to better distribute the pressure points



Product Reference	Description	Size	Distributor code
472333	Respireo™ N	S	KM213000
472334	Respireo™ N	M	KM213100
472335	Respireo™ N	L	KM213200
472336	Respireo™ N	XL	KM221700
472344	Respireo™ Primo N	S	KM223200
472345	Respireo™ Primo N	M	KM223300
472346	Respireo™ Primo N	L	KM223400
472347	Respireo™ Primo N	XL	KM223500

Spare parts:

472459	Respireo™ N/Primo N cushion x10	S	KM223600
472461	Respireo™ N/Primo N cushion x10	M	KM223700
472463	Respireo™ N/Primo N cushion x10	L	KM223800
472439	Respireo™ N/Primo N cushion x10	XL	KM221800
472412	Respireo™ Blue Non Vented elbow*	-	KM229109
472456	Generic Headgear	-	KM224000

Manufactured by

AIR LIQUIDE MEDICAL SYSTEMS S.r.l
Via Bisceglie, 66 - 20152 Milano (MI) - ITALY
PLANT: Via dei Prati, 62 - 25073 Bovezzo (BS) - ITALY
Tel. +39 030.201.59.11
Fax +39 030.209.83.29

Distributor Contact

AIR LIQUIDE MEDICAL SYSTEMS S.A.
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex
France

medicaldevice.airliquide.com

Hotline ALMS: +33 179 51 70 01

Air Liquide
HEALTHCARE

Respireo™ N / Primo N is a reusable nasal mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by a single adult patient (weighing more than 30 kg) who has been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

Medical devices class IIa CE0051. Please read the user manual.

* Medical device that changes the intended use of the mask into the following one:
Respireo™ Primo N with the Respireo™ Blue Non Vented elbow connector mounted is a reusable nasal mask without calibrated exhalation orifices, intended to be used at home or in a hospital setting by multiple adult patients (weighing over 30kg) in controlled non-invasive positive pressure ventilation (NIPPV).